



# SODIUM CITRATE TUBES FOR COAGULATION

Tubos con Citrato de Sodio para Coagulación / Tubes avec Citrate de Sodium pour Coagulation / Provette con Sodio Citrato per Coagulazione



## TEST TUBE WITH SODIUM CITRATE FOR COAGULATION

tubo con Citrato de Sodio para Coagulación / tube avec Citrate de Sodium pour Coagulation / provetta con Sodio Citrato per Coagulazione

## PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

## UNBREAKABLE

irrompible / incassable / infrangibile

## BLUE

azul / bleu / blu









## POLYETHYLENE

polietileno / polyéthylène / polietilene



## POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	22220	0,35 ml of Sodium x 3,15 ml of blood x ACL 2000-7000-9000	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32218	0,25 ml of Sodium x 2,25 ml of blood	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32247	0,3 ml of Sodium x 2,7 ml of blood	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32250	0,4 ml of Sodium x 3,6 ml of blood	13x75 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32210	0,25 ml of Sodium x 2,25 ml of blood	12x86 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi
	32254	0,5 ml of Sodium x 4,5 ml of blood	12x86 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi