

PETRI DISH WITHOUT VENTS

placa de petri sin vientos / boîte de petri sans ergots / piastra di petri senza ventilazione

MAXIMUM TRANSPARENCY

máxima transparencia / transparence maximale / massima trasparenza


PERFECT PLANARITY

planaridad perfecta / parfaite planéité / planarità perfetta



POLYSTYRENE

poliestireno / polystyrène / polistirolo

	CODE	TYPE	STERILITY	PACKAGING
	29140	ø 60 mm, H 14,2 mm	PRODUCED IN CLEANROOM ISO 6 producido en cleanroom ISO 6 / produit en cleanroom ISO 6 / prodotta in camera bianca ISO 6	BAG OF 10 PIECES bolsa de 10 piezas / en sachet de 10 pièces / in sacchetto da 10 pezzi
	29150	ø 60 mm, H 14,2 mm	STERILE estéril / stérile / sterile	BAG OF 10 PIECES bolsa de 10 piezas / en sachet de 10 pièces / in sacchetto da 10 pezzi