

Vacumed® Double wall Sodium Citrate test tube

FL Medical works to satisfy customers' needs, The company is continuously updating to improve the products and realize new projects.
At the basis of this item there are FL Medical values

Innovation
Quality
Efficiency

Main characteristics of the Vacumed® Double wall Sodium Citrate test tube are

INNER TUBE

Material: Polypropylene

Capacity: 1,8 ml or 2,7 ml

Polypropylene permits a better sodium citrate maintenance guaranteeing the correct additive/blood volume ratio

OUTER TUBE

Dimension: 13 x 75 mm

Material: PET

PET retains vacuum ensuring the blood volume/additive ratio

CAP


Dimension: 12 mm

Shield Material: Polypropylene

Stopper Material: Rubber

Colour: Light Blue

Main features of the cap are the perfect design and the high quality of materials which provide a perfect tightness



Code	Sodium Citrate	Aspiration Volume	Packaging
42450	3,2%	1,8 ml	100 pcs per rack 1.200 pcs per box
42451	3,2%	2,7 ml	100 pcs per rack 1.200 pcs per box
42490	3,8%	1,8 ml	100 pcs per rack 1.200 pcs per box
42491	3,8%	2,7 ml	100 pcs per rack 1.200 pcs per box