

Vacumed[®] tubes with Boric Acid

FL Medical works to satisfy customers' needs, the company is continuously updating to improve the products and realize new projects.
At the basis of this item there are FL Medical values

Innovation
Quality
Efficiency

The main characteristics of the Vacumed[®] tubes with boric acid are

Storage of the sample

Boric Acid allows to keep urine sample during warmest periods or when analysis is not possible in short time

Rapid Dissolution

Boric Acid tablet permits a rapid dissolution after sample collection

Results within 48 hours

Boric Acid allows to test the urine collected within 48 hours



Code	Tube dimensions	Aspiration Volume	Packaging
42972	13 x 75 mm	4 ml	100 pcs per rack 1.200 pcs per box
43976	13 x 100 mm	6 ml	100 pcs per rack 1.200 pcs per box
44970	16 x 100 mm	10 ml	100 pcs per rack 1.200 pcs per box