



TEST TUBES IN POLYSTYRENE FOR URINE

Tubos de ensayo en poliestireno para orina / Tubes en polystyrène pour urine /
Provette in polistirolo per urina



TEST TUBE CONICAL

tubo cónico / tube conique / provetta conica

WIDE MOUTH

embocadura ancha / large embouchure /
imboccatura larga

GRADUATED

graduado / gradué / graduata

MAXIMUM TRANSPARENCY

máxima transparencia / transparence maximale /
massima trasparenza

EXCELLENT RESISTANCE

excelente resistencia / excellente résistance /
eccellente resistenza



POLYSTYRENE

poliestireno / polystyrène / polistirolo

CODE	CAPACITY	TYPE	PACKAGING
25062	12 ml	17x105 mm	IN BULK a granel / en vrac / sfusa