



TEST TUBE CYLINDRICAL

tubo cilíndrico / tube cylindrique /
provetta cilindrica

EXCELLENT RESISTANCE

excelente resistencia / excellente résistance /
eccellente resistenza

STERILE

estéril / stérile / sterile

RED

rojo / rouge / rosso





POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène /
polipropilene

	CODE	CAPACITY	TYPE	PACKAGING
	25095	12 ml	16x100 mm (cylindrical)	INDIVIDUALLY WRAPPED embalado individualmente / sachet unitaire / in sacchetto singolo
	25096	12 ml	16x100 mm (conical)	INDIVIDUALLY WRAPPED embalado individualmente / sachet unitaire / in sacchetto singolo