



SEDI-RATE PIPETTE

pipeta Sedi-Rate / pipette Sedi-Rate / pipetta Sedi-Rate

GRADUATED FROM 0 TO 180 MM

graduada de 0 a 180 mm / gradué de 0 à 180 mm /
graduata da 0 a 180 mm

POLYSTYRENE

poliestireno / polystyrène / polistirolo

CODE	TYPE
50334	Sedi-Rate



TAKIVES IN GLASS PIPETTE WITH AUTOMATIC ASPIRATOR

pipeta Takives en vidrio con aspirador automático /
pipette Takives en verre avec aspirateur
automatique / pipetta Takives in vetro con
aspiratore automatico

GRADUATED FROM 0 TO 160 MM

graduada de 0 a 160 mm / gradué de 0 à 160 mm /
graduata da 0 a 160 mm

GLASS

vidrio / verre / vetro

CODE	TYPE
50315	Takives