



SODIUM CITRATE TUBES FOR ESR

Tubos con Citrato de Sodio para VSG / Tubes avec Citrate de Sodium pour VS /
Provette con Sodio Citrato per VES



SEDI-RATE TEST TUBE WITH 0,2 ML OF SODIUM CITRATE (3,8%) X 0,8 ML OF BLOOD (ESR)

tubo Sedi-Rate con 0,2 ml de Citrato de Sodio (3,8%) x 0,8 ml de sangre (VSG) / Sedi-Rate tube avec 0,2 ml de Citrate de Sodium (3,8 %) x 0,8 ml de sang (VS) / provetta Sedi-Rate con 0,2 ml di Sodio Citrato (3,8%) x 0,8 ml di sangue (VES)

PERFECT TRACEABILITY

trazabilidad perfecta / traçabilité parfaite /
tracciabilità perfetta

UNBREAKABLE

irrompible / incassable / infrangibile

PINK

rosa / rose / rosa



POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène /
polipropilene

CODE	CAPACITY	TYPE	PACKAGING
50335	0,8 ml	12x86 mm	RACK OF 50 PIECES Rack de 50 pièces / en portoir de 50 pièces / in rack da 50 pezzi