



**TEST TUBE WITH 0,1 ML OF SODIUM CITRATE (3,8%) X 0,9 ML OF BLOOD (COAGULATION)**

tubo con 0,1 ml de Citrato de Sodio (3,8%) x 0,9 ml de sangre (Coagulación) / tube avec 0,1 ml de Citrate de Sodium (3,8%) x 0,9 ml de sang (Coagulation) / provetta con 0,1 ml di Sodio Citrato (3,8%) x 0,9 ml di sangue (Coagulazione)

**PERFECT TRACEABILITY**

trazabilidad perfecta / traçabilité parfaite / tracciabilità perfetta

**UNBREAKABLE**

irrompibile / incassable / infrangibile



**POLYETHYLENE**

polietileno / polyéthylène / polietilene



**POLYPROPYLENE**

polipropileno / polypropylène / polipropilene

CODE	CAPACITY	TYPE	PACKAGING
22239	0,9 ml	12x56 mm	RACK OF 50 PIECES Rack de 50 piezas / en portoir de 50 pièces / in rack da 50 pezzi