



URINE CONTAINER WITH PAPER LABEL

frasco para orina con etiqueta de papel / pot pour urine avec étiquette en papier / contenitore per urina con etichetta cartacea

ASPIRATION SYSTEM FOR VACUUM TUBES

sistema de aspiración para tubos de vacío / système d'aspiration pour tubes sous vide / con sistema di aspirazione per provette sottovuoto

WITH SCREW CAP

con tapón roscado / avec bouchon à vis jointé / con tappo a vite







POLYETHYLENE

polietileno / polyéthylène / polietilene



POLYPROPYLENE

polipropileno / polypropylène / polipropilene

	CODE	CAPACITY	STERILITY	PACKAGING
	25014	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	IN BULK a granel / en vrac / sfusa
	25012	120 ml	PRODUCED IN CLEANROOM ISO 8 producido en cleanroom ISO 8 / produit en cleanroom ISO 8 / prodotto in camera bianca ISO 8	INDIVIDUALLY WRAPPED envuelto individualmente / sachet unitaire / in sacchetto singolo
	25016	120 ml	STERILE estéril / stérile / sterile	IN BULK a granel / en vrac / sfusa
	25015	120 ml	STERILE estéril / stérile / sterile	INDIVIDUALLY WRAPPED envuelto individualmente / sachet unitaire / in sacchetto singolo