




L-SHAPE SPREADERS

Asas a L / Étaleurs "L" / Aste a L



L-SHAPE SPREADER

asa a L / étaleur "L" / asta a L

STERILE

estéril / stérile / sterile



WHITE

blanca / blanc / bianca



POLYSTYRENE

poliestireno / polystyrène / polistirolo

	CODE	TYPE	PACKAGING
	27410	L-Shape	INDIVIDUALLY WRAPPED IN PEEL-PACK envuelto individualmente en peel-pack / sachet unitaire en peel-pack / in sacchetto singolo peel-pack
	27411	L-Shape	PACKET OF 5 PCS IN PEEL-PACK peel-pack de 5 piezas / sachet de 5 pièces en peel-pack / in sacchetto da 5 pezzi peel-pack